COD. 223678 Rev. 03 Data 2016/12 LINEA Carlevale

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Table with symbols and technical specifications: Lot number, Body diameter, Total diameter, Use before, Read instructions, Sterilized by steam, Dioptric power, Serial Number.

COD 223678 rev03 2016/12

REVIEW PREMIUM IOL

FIL SSF Carlevale lens

Sutureless Scleral Fixation

1. DESCRIPTION

Carlevale hydrophilic Intraocular lenses are manufactured from PolyHema 25% H2O with ultraviolet absorber, are available in several models, for the implantation in the posterior chamber according to standard surgical procedures according to the traditional surgical methods of scleral fixation iol's implant.

2. PACKAGING


Carlevale hydrophilic Intraocular lenses are stored in bidistilled water free from pirogens, steam sterilised in autoclave and placed in vials. Each sterile lens is enclosed in its own lens "Holder Folder" inside a vial and a double sterilisation pouch. The external surfaces of the outer pouch are not sterile. The pouch including the product, the sticky labels, the Information form and the patient card are inclosed in a single package. The package includes a patient card that must be given to the patient. The patient must keep this card and show it to his or her doctor when undergoing routine checkups.

3. CHARACTERISTICS

Carlevale hydrophilic Intraocular lenses are available in a standard diopter range from +10.0 diopters to +30.0 diopters, with a 0.50 diopter increment from +15.0 diopters to +25.0 diopters, and a 1.00 diopter increment for all other values. Upon request, some models are available from -10.0 diopters to +10.0 diopters. The resolution efficiency of these lenses is 60% or greater Soleko ultraviolet absorbing. Intraocular lenses have an ultraviolet and visible light transmittance of 10% at 374 nm for a +10.0 diopter lens and 10% at 386 nm for a +30.0 diopter lens (see Ultraviolet Transmittance).

4. INDICATIONS

The implant of hydrophilic Carlevale iol is suitable for the correction of aphakic eyes in which there is not capsular bag/ support. Carlevale iol's are designed to be implanted in the sulcus without capsular bag by anchoring of special self blocking plugs on the bed of scleral pocket and covered by scleral flaps (minimum suggested dimensions 3.5 mm X 3.5 mm).

5. CONTRAINDICATIONS

Patients with one of the following particularly cannot be suitable for the implant of Carlevale IOL.

- 1. Inflammation of anterior and posterior segment.
2. Recurrent anterior or posterior segment inflammation of unknown aetiology.
3. Patients on whom the Intraocular lens may interfere with the ability to observe, diagnose or treat posterior segment disease.
4. Surgical difficulties at the time of cataract extraction that might increase the potential for complications (e.g.: persistent bleeding, significant iris damage, uncontrollable positive pressure, or damage due to significant vitreous prolapsis or loss).
5. Glaucoma with little pharmacological control, neovascular glaucoma, glaucoma with through draft and / or vascular drainage systems.
6. Corneal endothelial dystrophy.
7. Proliferative diabetic retinopathy.
8. Severe abnormalities and clotting dysfunction.
9. Abnormalities, trauma and thinning of the sclera.

6. COMPLICATIONS AND ADVERSE EFFECTS

As in any surgery, risks are present, much more with this type of iol that is implanted in case patients have already had complications of rupture or capsular removal: dislocation or iol tilting, tearing of self blocking plug and/or haptics' parts of iol both during loading of iol in the cartridge and during pulling out plug from eyeball on the bed of scleral flap, hyphema and / or emovitreo, acute or chronic inflammation, endothelial injury, endophthalmitis, extrusion of self blocking plugs from the sclera and conjunctiva, chronic trauma of the root and/or rear portion of the iris, pigment dispersion, secondary glaucoma, transitional overtone, retinal detachment, vitritis, cystoid edema, pupillary membrane, iris prolapse, ipopion. Adverse reactions (hypopyon, Intraocular infection, acute corneal decompensation and/or secondary surgical intervention) and/or potentially sight-threatening complications that may reasonably be regarded as lens related and that were previously expected in nature, severity or degree of incidence should be promptly reported to Soleko.



7. WARNINGS

The security of Carlevale scleral fixation iol is not tested in patients with very high myopia and/or diseases, trauma or thinning of scleral limbus, tube draining systems, neovascular glaucoma, glaucoma with through draft and / or vascular drainage systems and glaucoma with control pharmacology precariously.

OTHER WARNINGS: The Carlevale iol is designed and projected to be implanted exclusively in posterior chamber making two sclera flaps on limbus basis with minimum dimension of 3.5 x 3.5 mm on the axis 0-180° and self blocking plugs have to be grasped with two 23/25° G vitrectomy forceps angulated with crocodile tip. The 23/25 G angulated forceps grasp the self blocking plugs and pull/take them out from the bed of sclera flap through an hole made by 23G sclerotome far 1.3 mm from limbus. As soon as the self blocking flaps are placed on the sclera bed, they have to be covered by sclera flaps, that together conjunctiva have to be sutured.

HYDROPHILIC

HYDROPHILIC ACRYLIC INTRAOCULAR LENSES • INFORMATION SHEET

8. CALCULATION OF LENS POWER

The physician should determine the power of the lens to be implanted. Calculations should estimate from the refractive error or calculations based on the corneal radius, depth of the anterior chamber and the paraxial length of the eye according to the several formulas published or existing in literature. The stated Constant, considering the various factors joining in the calculation, is estimated it is similar to an IOL for capsular bag.

9. PRECAUTIONS

The Carlevale intraocular lens must not be resterilised. Opening the sterile pouch requires an immediate use of the Carlevale intraocular lens. Check the shelf-pack integrity prior to use. Any damaged shelf-pack should be returned to Soleko with its contents, in accordance with the return lens policy described below. The lenses must be kept at a temperature between 20° and 25°C. If kept at lower temperature, they must be reconditioned, for at least 3 hours, at a temperature between 20 and 25°C.

- Warning: check the surface of the lens by looking through the vial. In case of opacity of the surface, do not use the lens and recondition it as above specified.
• In case the opacity is noted when the vial has already been opened, do not use the lens.

10. LENS RETURN POLICY

Please, contact Soleko for any returns: Tel. 0039 02 37924126.

11. INSTRUCTIONS FOR USE

There are various surgical procedures that can be used for the Intraocular lens implantation. Following inspection, the lens will be ready for insertion. During the surgical operation, it may be necessary to use a viscoelastic or other transitional medium for implantation.

- 1. Before opening the package, check the label for model, dioptric power and expiration date.
2. Inspect the external pouch to insure that sterility has been maintained.
3. Open the sterile pouch in a sterile field and remove the vial. Again confirm the lens dioptric power and the stated model.
4. Open the foil and extract the lens holder from the vial.
5. Grab the lens with Lens Inserting tweezers.
6. Examine the lens placed in the sterile forceps. Once more verify the stated model. Hydrate the lens using a sterile Intraocular irrigation solution (Balanced Saline Solution).
7. Hold the lens in the Lens Inserting tweezers. Either insert the lens in the patient eye or transfer the IOL into a disposable cartridge. During the loading in the cartridge, place the carlevale lens considering the indicator of way (showed by an arrow on the design) so that the iol will be in the right way after injecting it into the eyeball and catching it through a 23/25 G curved crocodile vitrectomy forceps.
8. Inspect the optic and haptic for defects and scratches that may have occurred during transportation and/or manipulation. As the lens and packaging materials may pick up an electrostatic charge, the lens should be carefully examined to ensure that no particles are attached.

EXPIRATION DATE

The expiration date refers to the sterility of the packaging. Do not implant the lens after this date.

Table with symbols and technical specifications: Lot number, Body diameter, Total diameter, Use before, Read instructions, Sterilized by steam, Dioptric power, Serial Number.