

Capsular Tension Rings

Flexible Insertion, Stable Result

Highly flexible due to the compression molding technology. These PMMA Capsular Tension Rings are indicated for the stabilization of weakened, broken or missing zonules preventing IOL decentration after capsular shrinking.



Capsular Tension Ring | CTR Model 275 12/10



Compression: From 12 mm to 10 mm PMMA | Flexible Ring Material: Clear PMMA Available in: Overall Ø: 12 mm

Capsular Tension Ring | CTR Model 276 13/11



Compression: From 13 mm to 11 mm Material: PMMA | Flexible Ring Available in: Clear PMMA Overall Ø: 13 mm

Features & Benefits

- Circular expansion and stabilization of the capsular bag
- Safe IOL centration in eyes with zonular dehiscence
- Prevents IOL decentration after capsular shrinking
- Stabilized conditions during (phaco emulsification) surgery
- · Reduced risk of capsular fibrosis



Recommended to be used in combination with the ON 385 EASY Control™ Micro Inserter











